

Considerations on Design Parameters


Preliminary Considerations:

An integrating sphere collects light through its input port and reflects this light from its diffuse, but highly reflective surface.

In this manner, the beam is spatially homogenized and it loses all polarization properties. This makes the integrating sphere an ideal instrument for preparation of a laser beam for power measurement.

In addition, the total power exiting a given port is proportional to the size of that port. In this manner, the integrating sphere can be used to attenuate the power of a laser beam making measurement with a photodiode accessible.

Since photodiodes are high speed devices, an integrating sphere is a much faster power measurement device than a thermopile.



The power reaching the photodiode mounted in an integrating sphere is given by:

$$P_{pd} = P_{laser} \times \frac{A_{pd}}{A_{sphere}} \times \frac{\rho}{1 - \rho(1 - a)}$$

Whereby:

P_{pd} = power onto the photodiode

A_{sphere} = surface area of the sphere

P_{laser} = power entering the sphere

ρ = reflectivity of the sphere surface


A_{pd} = surface area of the photodiode

$a = \frac{\text{combined surface area of all ports}}{\text{surface area of the sphere}}$

From this we see that a small value of „a“ increases the sensitivity of the sphere. It also has the effect of improving the homogenizing property of the sphere leading to more reproducible measurements.

The ratio $\frac{A_{pd}}{A_{sphere}}$ allows us to adjust the sensitivity of the measurement system to match a particular application.

Light may be coupled from the integrating sphere to an external photodiode via an optical fibre. This further attenuates the power allowing measurement of much higher laser powers. Furthermore, since the photodiode is not mounted in the integrating sphere, it will not heat up at high power. This conserves the integrity of the calibration of the system. Artifex Engineering calibrates such systems as a complete chain to ensure accuracy of the measurement system.



Standard Options (Polymer Spheres):

Spheres with integrated photodiode and SMA fibre port:

Polymer Integrating Spheres P10



Polymer Integrating Spheres P20



- High reflectivity from 250-2500nm
- Very homogeneous
- Power handling up to 100W

Polymer Integrating Spheres SP20



Polymer Integrating Spheres SP50



Polymer Integrating Spheres SP100




Standard Options (Gold Spheres):

- High reflectivity from 700-12000nm
- Power handling up to 12kW


25mm Gold Integrating Spheres
G25



65mm Gold Integrating Spheres
G65





Liquid Cooled, 100mm
Gold Integrating Spheres
G100L



Liquid Cooled, 200mm
Gold Integrating Spheres
G200L




Standard Options (Detector Ports):



SP- Series with rotatable Azimuth


Standard Options (Photodiode Ports):



Choices of Photodiodes

- UVS (Si): 250-1100nm
250-950nm preferred
- VIGA (InGaAs): 400-1600nm
- x2.0IGA (x-InGaAs): 800-2000nm
- x2.2IGA (x-InGaAs): 1000-2200nm
- x2.5IGA (x-InGaAs): 1300-2500nm

Standard Options (Entrance Ports):



- Proprietary design for very low back reflection. This means the adapter will not influence the calibration.
- FC and SMA available.

Take Aways:

- Integrating spheres fitted with photodiodes are accurate power measurement devices.
- Integrating spheres are much faster than thermopiles.
- Spheres can be configured for almost any application.
- Power measurement over 8 decades possible with a single device.
- Power range: 1 μ W (min.) up to 12kW (max.).
- Wavelength range (calibrated):
250-2490nm (max. 100W)
or
700-2490nm (max. 12kW).